



Specifications	Value	Specifications	Value
X-ray Generator		Mobile Stand	
CubeX 28 HF X-ray Generator	2.7 kW, 40 - 100 kVp, 40 mA Max 28 lbs (13 kg)	Weight	92 lbs (41.7 kg) with generator and cables
Controls	Integrated controls with ExamVue Duo Software Inverted Control panel for use on a stand or table	Maximum SID to Floor	79" (200 cm)
Panel		Large Wheels	10" (25.4 cm)
ExamVue Duo EV14W	14" x 17" Wireless Detector Assembled in the US 140 µm pixel pitch, 3072 x 2560 pixels active area Direct Deposit CSI with AED Dual Band Wireless Support (2.4 Ghz and 5 Ghz) 3.4 kg (7.5 lbs) Quick connect magnetic cable Includes 2 batteries and battery charger	Folded Dimensions	27" (61 cm) x 57" (137 cm) x 27" (61 cm)
Panel Protection Options		Integrated Cabling	Exposure Switch Extension
Mobile DR Panel Holder	14" x 17" with Single Handle Long Edge	Acquisition Software	
Full Panel Encasement With Snap-On Grid	103 LPI 14" x 17" with Aluminum Frame for Mobile	ExamVue Duo	Acquisition Software for Medicine
		Acquisition Laptop	* Specs subject to change without notice Windows 10 PRO 16 GB RAM, 512 GB SSD Hard Drive i5 10th GEN Processor 15.6" display
		Support	
		Remote Service Support Software	Team Viewer Remote Access Software Included (allows for WEB-Based Remote dial-up support)

Contact us:

3555 Veterans Memorial Hwy, Unit D
Ronkonkoma, NY 11779
sales@jpihealthcare.com
(516) 513-1330 Ext. 3



www.jpihealthcare.com