



WORLD'S BEST SELLING BATTERY GENERATOR VET-20BT

- ➔ 15 LBS WEIGHT
- ➔ BATTERY OPERATED
- ➔ 90kV, 20mA, 1.6kW OUTPUT



Max. Output(kW)	1.6
kV Range per 1kV step	50-90
mA Range/No. of step	8, 10, 15, 20
mAs Range/No. of step	0.4-20/25
Tube Focal Spot(mm)	1.2 x 1.2
Target Angle	19'
Dual Laser Pointer	Standard
Battery Charger	16.8V/1.5A
Charging Time	3.5 hours
No. of APR	15
Packing Dimension	51 x 36 x 37 (37lbs)